Q CREMA EDTA GEL

(E)

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(E N)

ROOT CANAL CONDITIONER

Indicaciones/properties: • The ethylene diamine tetra acetic acid and its derivatives have a strong chelating power over various ions. This property gives them many applications. In Odontology they are widely used as dentin softeners to make easier the dentin remove.

• In core preparations Q Crema remove easily the smear layer. • In canal instrumentations placing Q Crema in the canal helped by a K file or similar, the dentin gets softer making easier and faster the procedure. • Q Crema contains Carbamide peroxide, which produces bubbles when the canal is irrigated with hypochlorite in the usual way. This effect helps cleaning out all rests of the instrumentation such as debris, smear layer, tissues, etc. • Without any fatty or non water soluble compounds a single rinse avoid leaving any rest inside the canal.

Instructions for use: • Open the access to the pulp chamber in the usual way. • Place Q Crema on a file and introduce it inside the canal, its recommended to get into the canal step by step, avoiding compressing the results

of the instrumentation towards the apex and trough it

- Only 1/3 of the canal length should be acceded in the first step. • Then irrigate with hypochlorite; afterwards more length can be reached. • The bubbles produced by the contact between the hypochlorite and the Q Crema make the debris, smear layer and all rests gets out from the canal.

Precautions: • A rubber dam is recommended to isolate the area. • Do not leave rest of Q Crema inside the canal and avoid the product pass through the apex as it can produce discomfort. Even if it happens no toxic effects should appear. • After the instrumentation the canal must be rinsed and dried completely with paper points. • Caution to reactions allergic to EDTA, or to carbamide peroxide. • Keep the product tightly closed.

Composition: • EDTA (Na₂) 15% • Carbamide peroxide 10% • Excipients.

Storage: • The product must be stored refrigerated between 4-18°C.

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