EN **Instructions for use**
Die silicone

Product description:
This **die silicone** is an addition-curing silicone for the extraoral fabrication of dental models. Indirect composite restorations can be fabricated chairside on this silicone model.
The **die silicone's** special curing formula allows you ample time to apply the material. The material is fully cured after 4 minutes, at which time you can then start modelling the restorations.

Indications:
Extraoral fabrication of cavity models and models of individual, or groups of teeth

Method of application:

Preparation:
Insert the cartridge into a type 2 **VOCO** dispenser or comparable dispenser. Due to technical reasons there may be slight differences in the fill levels of the two cartridge chambers prior to initial use.

In such a case, remove the cartridge cap and carefully extract material until it is delivered evenly from both openings. Then attach a type 20 mixing tip and lock with a 90° turn clockwise. The material is extruded by pressing the lever and is automatically mixed in the correct ratio.

The attached mixing tip is intended for single use only. After use, store the cartridge tightly sealed. Check that material flows freely through the openings before reusing.

Processing:
To fabricate a cavity model, or a model of an individual tooth or several teeth, first take an impression of the prepared area with an alginate. Silicones or polyethers cannot be used for this purpose. Clean and dry the impression in the usual way. Ensure all drops of moisture have been removed as these will lead to artefacts in the model and thus in the restoration. Then fill the impression with the **die silicone** within the working time.

Before use, discard the material (roughly the size of a bean) that is initially extruded from the mixing tip.

Starting at the deepest point, place the material directly into the impression and fill the cavity from the bottom upwards leaving the tip immersed. When applying the material, ensure that it flows evenly to avoid air bubbles.

Working time at room temperature (23 °C) is at least 30 seconds. Higher temperatures reduce and lower temperatures extend the working and setting times. The completed silicone model can be removed from the impression 4 minutes after commencement of mixing.

After the material has set it is extremely hard, enabling you to start modelling the composite restoration immediately. Follow the instructions for use of the composite regarding preparation, processing and curing.

Indirect composite restorations can then be adhesively attached in the usual way, following the respective manufacturer's instructions for use.

Information, precautionary measures:

- This **die silicone** contains a platinum catalyst. The **die silicone** should therefore not be used in case of known hypersensitivities (allergies) to this ingredient.
- Do not use in combination with condensation-curing silicones, polyethers or materials containing eugenol.
- Setting may be impaired by the use of latex gloves, hand creams, cleaning agents, mouth rinses, retraction solutions, materials containing methacrylates, or residues of other materials. These should therefore be removed before applying the material.
- Avoid incompletely cured silicone material coming into contact with clothing.
- Please note that any residues of material left in the mixing tips will cure and set.
- The warranted product characteristics can only be achieved when using the original Sulzer Mixpac mixing tips as approved by **VOCO**. All the mixing tips for use with our respective products are available from **VOCO**.

Storage:
Store at temperatures of 4 °C to 23 °C. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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