

mectron



→ MECTRON  
ULTRASONIC INSERTS




# → ULTRASONIC INSERTS

→ MULTIPIEZO

## → A LARGE VARIETY OF APPLICATIONS

Scaling, perio, endo, restorative – Mectron offers a wide range of applications. Mectron's high-quality technique guarantees optimal treatment power and utmost security in each single unit.

This remarkable quality is the result of the successful combination of fine electronic concept and a special tips design.



→ COMBI touch



→ MICROPIEZO



→ COMPACT PIEZO LED



→ COMPACT PIEZO P2K



# → THE INSERTS




Mectron offers an unequalled variety of inserts – for every indication, every treatment, and every situation you will face in your daily practice. More than 30 different inserts are available.

# → TORQUE WRENCH K10



## → TORQUE WRENCH K10 – PROTECTS YOUR HANDPIECE AND INSERTS

- sterilisable
- inserts are supplied with their own dynamometric wrench





# CONTENTS

---

→ 8-14 INSERT SETS

8/9 Set maintenance/Set scaling

10/11 Set perio maintenance/Set perio treatment

12/13 Set implant cleaning/Set ortho clean-up

14 Set endo revision

→ 16-22 SCALING INSERTS

16/17 Insert S1/Insert S1-S

18/19 Insert S2/Insert S3

20/21 Insert S4/Insert S5

22 Insert S6

→ 24-27 PERIO UNIVERSAL INSERTS

24/25 Insert P1/Insert P2

26/27 Insert P3/Insert P4

→ 28-35 PERIO ANATOMIC INSERTS

28/29 Insert P10/Insert P11

30/31 Insert P12/Insert P13

32/33 Insert P14/Insert P15

34/35 Insert P16R/Insert P16L

→ 36-37 IMPLANT CLEANING INSERT

→ 38-39 ORTHODONTIC CLEAN-UP INSERTS

38/39 Insert D05120/Insert D07102 + tip D06007

→ 40-45 ENDO REVISION INSERTS

41 Insert ER1

42/43 Insert ER2/Insert ER3

44/45 Insert ER4/Insert ER5

→ 46-49 ENDO RETRO INSERTS

47 Insert RD3

48/49 Insert RD4/Insert RD5

→ 50-53 EXTRACTION INSERTS

51 Insert ME1

52/53 Insert ME2/Insert ME3

→ 54-55 CROWN PREP TIPS

→ 56-59 APPROPRIATE INSERT SETTINGS

→ 60-61 REFERENCE NUMBERS/TIP CARD

# → SET MAINTENANCE

→ EQUIPPED WITH:

- 1 insert S1
- 1 insert S2
- 1 insert S5
- 1 insert P2
- 1 insert P3
- 5 torque wrenches K10





# → SET SCALING

## → EQUIPPED WITH:

- 1 insert S1
- 1 insert S2
- 1 insert S3
- 1 insert S4
- 1 insert S5
- 5 torque wrenches K10



# → SET PERIO MAINTENANCE

## → EQUIPPED WITH:

- 1 insert P10
- 1 insert P11
- 1 insert P12
- 1 insert P13
- 1 insert P14
- 5 torque wrenches K10



# → SET PERIO TREATMENT



- EQUIPPED WITH:
- 1 insert P2
  - 1 insert P3
  - 1 insert P15
  - 1 insert P16R
  - 1 insert P16L
  - 5 torque wrenches K10



INSERT SETS


# → SET IMPLANT CLEANING

→ EQUIPPED WITH:

- 1 insert ICS
- 1 torque wrench K10
- 5 Implant cleaning tips IC1



# → SET ORTHO CLEAN-UP\*



- EQUIPPED WITH:
- 1 insert D05120
  - 1 insert D07102
  - 2 torque wrenches K10
  - 5 D06007 tips



\* developed in collaboration with Prof. Gianna Maria Nardi


# → SET ENDO REVISION

→ EQUIPPED WITH:


- 1 insert ER1
- 1 insert ER2
- 1 insert ER3
- 1 insert ER4
- 1 insert ER5
- 5 torque wrenches K10









- IDENTIFICATION  
scaling
- MORPHOLOGY  
universal curette with semicircular diameter
- SURFACE  
brightly polished medical grade stainless steel
- TREATMENT  
for considerable tartar removal









- **IDENTIFICATION**  
perio/scaling
- **MORPHOLOGY**  
long universal curette with triangled slightly curved surface
- **SURFACE**  
brightly polished medical grade stainless steel
- **TREATMENT**  
for supragingival tartar removal and subgingival concretions and biofilm elimination








- **IDENTIFICATION**  
scaling
- **MORPHOLOGY**  
universal curette with triangled, slightly curved surface
- **SURFACE**  
brightly polished medical grade stainless steel
- **TREATMENT**  
efficient in the interdental spaces and posterior surfaces





- IDENTIFICATION  
scaling
- MORPHOLOGY  
flat, with rounded edges
- SURFACE  
brightly polished medical grade stainless steel
- TREATMENT  
for considerable supragingival  
tartar removal





→ IDENTIFICATION

scaling

→ MORPHOLOGY

universal curette with 45° angled, triangled and slightly curved surface

→ SURFACE

brightly polished medical grade stainless steel



→ TREATMENT


for tartar removal on mesial and distal interdental surfaces in the posterior area






- **IDENTIFICATION**  
scaling
- **MORPHOLOGY**  
similar shape to S1 insert  
but longer and thinner
- **SURFACE**  
brightly polished medical grade stainless steel
- **TREATMENT**  
for gentle supra- and subgingival  
tartar removal and for gingivitis






- **IDENTIFICATION**  
scaling
- **MORPHOLOGY**  
contra-angled universal curette with  
semicircular diameter
- **SURFACE**  
brightly polished medical grade stainless steel
- **TREATMENT**  
powerful insert (twice the power of S1 insert)  
for considerable tartar removal








- IDENTIFICATION  
perio
- MORPHOLOGY  
long, straight working tip with a circular surface
- SURFACE  
brightly polished medical grade stainless steel
- TREATMENT  
for concrements and biofilm  
removal from root surfaces








- IDENTIFICATION  
perio
- MORPHOLOGY  
curved working tip with a circular surface
- SURFACE  
brightly polished medical grade stainless steel
- TREATMENT  
for concrements and biofilm  
removal from root surfaces





→ IDENTIFICATION

perio

→ MORPHOLOGY


slim and smooth like a probe

→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

supra and subgingival concretions and  
biofilm removal from interproximal spaces





→ IDENTIFICATION

perio

→ MORPHOLOGY



extra long, straight working tip  
with a circular surface


→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for concrements and biofilm removal  
from root surfaces





→ IDENTIFICATION

perio

→ MORPHOLOGY

extra long, curved working tip  
with a circular surface

→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for concrements and biofilm removal  
from root surfaces





→ IDENTIFICATION

perio

→ MORPHOLOGY


right (inclination of 15°) curved insert  
with round tip


→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for gentle subgingival concrements removal





→ IDENTIFICATION

perio

→ MORPHOLOGY

left (inclination of 15°) curved insert  
with round tip

→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for gentle subgingival concretions removal





→ IDENTIFICATION

perio

→ MORPHOLOGY


right (inclination of 15°) curved insert with spherical (Ø = 0.8 mm) round tip


→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for concrements and biofilm removal in furcations and concavities





→ IDENTIFICATION

perio

→ MORPHOLOGY

left (inclination of 15°) curved insert with spherical ( $\varnothing = 0.8 \text{ mm}$ ) round tip

→ SURFACE

brightly polished medical grade stainless steel

→ TREATMENT

for concrements and biofilm removal in furcations and concavities







→ IDENTIFICATION

perio

→ MORPHOLOGY



universal curette shape with round surface

→ SURFACE


brightly polished medical grade stainless steel

→ TREATMENT

supra- and subgingival concretions and biofilm removal. Recommended for the treatment of canine and anterior teeth. Debridement of deep periodontal pockets.



# → P16R



→ IDENTIFICATION

perio

→ MORPHOLOGY


right angled periodontal curette

→ SURFACE

brightly polished medical grade stainless steel



→ TREATMENT

subgingival concretions and biofilm removal from furcations and deep pockets. Recommended for the treatment of molars and premolars.






- IDENTIFICATION  
perio
- MORPHOLOGY  
left angled periodontal curette
- SURFACE  
brightly polished medical grade stainless steel
- TREATMENT  
subgingival concretions and biofilm removal from furcations and deep pockets. Recommended for the treatment of molars and premolars.



# → ICS + IC1

## → INSERT FOR IMPLANT CLEANING – MAXIMAL TREATMENT EFFICIENCY WITH THE HIGHEST PATIENT COMFORT

- tip IC1 – 100% PEEK (high tech plastics) – for the maximum respect of the implant surface
- the special shape of the tipholder ICS allows an optimal access and a perfect view even in posterior regions
- tipholder ICS and tip IC1 are both sterilisable and reusable



### → IMPLANT CLEANING SET S

The set consists  
of a tipholder ICS,  
1 torque wrench K10  
and 5 tips IC1



→ IDENTIFICATION

perio

→ MORPHOLOGY



102° angled tipholder with conical PEEK tip (working length 10 mm, Ø top 0.8 mm, Ø base 2 mm)

→ SURFACE

brightly polished medical grade stainless steel for tipholder ICS and 100% PEEK for IC1

→ TREATMENT

efficient and gentle cleaning of implants and restorations




→ ICS



→ IC1

implant cleaning tip IC1 can be easily tightened without special tools and is fully autoclavable

# → D05120



## → IDENTIFICATION

restorative

## → MORPHOLOGY

angled insert with  $\varnothing 0,60 > \varnothing 1,50$ , working length 9.5 mm, D30 diamond coated

## → SURFACE

brightly polished medical grade stainless steel D30 diamond coated

## → TREATMENT

safe and effective removal of the residual adhesive post orthodontic therapy (fixed, removable/invisible)





# D07102 + D06007



## → IDENTIFICATION

restorative

## → MORPHOLOGY


102° angled tipholder with conical PEEK tip  
(working length 11.5 mm, Ø top 0.92 mm,  
Ø base 2.4 mm)

## → SURFACE

brightly polished medical grade stainless steel for  
tipholder D07102 and PEEK for tip D06007

## → TREATMENT

gentle polishing post removal of residual adhesive



→ D07102



→ D06007

The tip D06007 made of PEEK is gentle on the enamel surface

## → ENDO REVISION WITH ER INSERTS


Exposure and removal of root pins, fast and effective removal of calcification in the coronal third of the root, removal of fractured root inserts – the ER inserts cover the complete spectrum of endodontic revision treatment.








- IDENTIFICATION  
endo revision
- MORPHOLOGY  
angled shaft with small, spherical (Ø 1.7 mm)  
and D30 diamond coated tip
- SURFACE  
titanium nitride coat
- TREATMENT  
exposure of root canal pins





→ IDENTIFICATION

endo revision

→ MORPHOLOGY

angled insert with 0.7 mm Ø, working length up to 10 mm, the last 5 mm D30 diamond coated

→ SURFACE

titanium nitride coat

→ TREATMENT

location of concealed or calcified root canal entrances, removal of restoration materials, calcification and fractured inserts in the coronal third of the root canal





→ IDENTIFICATION

endo revision

→ MORPHOLOGY


angled insert with 0.6 mm  $\varnothing$ , working length up to 10 mm, the last 5 mm D30 diamond coated

→ SURFACE


titanium nitride coat

→ TREATMENT

presentation of calcified root canals entrances, removal of root elements and fractured inserts in the coronal third of the root canal



# ER4



## → IDENTIFICATION

endo revision

## → MORPHOLOGY

angled insert with 0.6 mm  $\varnothing$ , working length up to 20 mm, the last 5 mm D30 diamond coated


## → SURFACE

titanium nitride coat


## → TREATMENT

presentation of calcified root canal entrances, removal of root elements and fractured inserts in the coronal third of the root canal






- IDENTIFICATION  
endo revision
- MORPHOLOGY  
angled insert with 0.5 mm  $\varnothing$ ,  
working length up to 24 mm
- SURFACE  
titanium nitride coat
- TREATMENT  
removal of fractured inserts in the medial  
and apical third of the root canal




## → ENDO RETRO INSERTS – MINIMALLY INVASIVE RETROGRADE ROOT CANAL PREPARATION


Endo retro inserts assure a satisfying solution to the root canal access issue. Their thin and 90° angled structure allows a conservative cut through a minimum wide bone window.

These advantages facilitate wound healing and spare the patient useless pain.






- IDENTIFICATION  
endo retro
- MORPHOLOGY  
angled shaft, with 90° angled working tip  
and fine diamond coating (D30)
- SURFACE  
titanium nitride coat
- TREATMENT  
root canal cleaning in the frontal area



# RD4



→ IDENTIFICATION

endo retro

→ MORPHOLOGY


angled shaft, with 90° right angled working tip and fine diamond coating (D30)

→ SURFACE


titanium nitride coat

→ TREATMENT


root canal cleaning in the molar area







- IDENTIFICATION  
endo retro
- MORPHOLOGY  
angled shaft, with 90° left angled working tip  
and fine diamond coating (D30)
- SURFACE  
titanium nitride coat
- TREATMENT  
root canal cleaning in the molar area



## → EXTRACTION INSERTS – EXCLUSIVE FOR UNITS WITH PULSE MODE


Even if the primary objective of dentistry is the preservation of natural teeth, sometimes extractions cannot be avoided. In these cases, maximum safety is the most important thing during treatment.

The specially shaped inserts ME1, ME2 and ME3 – with a blade thickness of only 0.35 mm – perfectly fulfill the anatomical demands of non surgical extractions. They offer maximum intra-operative sensitivity, allowing a great preservation of alveolar bone and periodontal tissues. Furthermore, the physiological solution guarantees maximum visibility thanks to the cavitation effect which leads to a temporary hemostasis of blood vessels.


→ EXCLUSIVE FOR UNITS WITH PULSE MODE




- IDENTIFICATION  
extractions
- MORPHOLOGY  
double-edged blade with a thickness of only 0.35 mm
- SURFACE  
titanium nitride coat
- TREATMENT  
atraumatic teeth removal



EXTRACTION INSERTS




- IDENTIFICATION  
extractions
- MORPHOLOGY  
left angled, double-edged blade  
with a thickness of only 0.35 mm
- SURFACE  
titanium nitride coat
- TREATMENT  
atraumatic teeth removal




→ EXCLUSIVE FOR UNITS WITH PULSE MODE

→ ME3\*



- IDENTIFICATION  
extractions
- MORPHOLOGY  
right angled, double-edged blade  
with a thickness of only 0.35 mm
- SURFACE  
titanium nitride coat
- TREATMENT  
atraumatic teeth removal



EXTRACTION INSERTS

# CROWN PREP TIPS

## CROWN PREP TIPS – ATRAUMATIC AND PRECISE

The crown prep tips are dedicated to the preparation and finishing of subgingival margins. Due to the particular morphology, the tip could be used in contact with the margin without damaging the surrounding soft tissues. The preparation margin could be positioned subgingivally giving more comfort to the patient.


The special shape of the tipholder allows an optimal view of the preparation field. The tipholder brings the crown prep tips to an elliptic movement and enables therefore a circular preparation of the tooth.




→ dynamometric wrench K7



→ key AB1



→ tip holder DB1 with crown prep tips



→ crown prep tip TF12D60



→ crown prep tip TA14D60

→ VERTICAL PREPARATION TECHNIQUE

→ LENGTH 10 MM → DIAMOND COATING

| diameter | D90     | D60     |
|----------|---------|---------|
| 1.2 mm   | TF12D90 | TF12D60 |
| 1.6 mm   | TF16D90 | TF16D60 |

→ CHAMFER TECHNIQUE

→ LENGTH 10 MM → DIAMOND COATING

| diameter | D120     | D90     | D60     |
|----------|----------|---------|---------|
| 1.2 mm   |          | TA12D90 | TA12D60 |
| 1.4 mm   | TA14D120 | TA14D90 | TA14D60 |
| 1.6 mm   | TA16D120 | TA16D90 | TA16D60 |

# →→→ APPROPRIATE INSERT SETTINGS



| function | endo perio/scaler restorative | endo perio/scaler restorative | endo perio/scaler prosthetics |  |
|----------|-------------------------------|-------------------------------|-------------------------------|--|
|----------|-------------------------------|-------------------------------|-------------------------------|--|

## SCALING

|                        |       |       |       |        |
|------------------------|-------|-------|-------|--------|
| S1-S                   | 1-2** | 1-2** | 1-2** | 0-30%  |
| S1, S2, S3, S4, S5, S6 | 1-6** | 1-5** | 1-3** | 0-100% |

## PERIO

|   |       |       |       |       |
|---|-------|-------|-------|-------|
| P1, P2, P3, P4, P10, P11, P12, P13, P14 | 1-3** | 1-3** | 1-3** | 0-50% |
| P15, P16R, P16L                         | 1-3** | 1-3** | 1-3** | 0-50% |

## IMPLANT CLEANING

|           |       |       |       |       |
|-----------|-------|-------|-------|-------|
| ICS + IC1 | 1-2** | 1-2** | 1-2** | 0-30% |
|-----------|-------|-------|-------|-------|

## ORTHO CLEAN-UP

|                         |     |     |   |         |
|-------------------------|-----|-----|---|---------|
| D05120, D07102 + D06007 | 3-6 | 3-5 | 3 | 50-100% |
|-------------------------|-----|-----|---|---------|





| function             | endo perio /scaler restorative | endo perio /scaler restorative | endo perio /scaler prosthetics |         |
|----------------------|--------------------------------|--------------------------------|--------------------------------|---------|
| <b>ENDO REVISION</b> |                                |                                |                                |         |
| ER1                  | 1-6                            | 1-5                            | 1-3                            | 0-100%  |
| ER2, ER3, ER4, ER5   | 1-3                            | 1-3                            | 1-3                            | 0-50%   |
| <b>ENDO RETRO</b>    |                                |                                |                                |         |
| RD3, RD4, RD5        | 1-3                            | 1-3                            | 1-3                            | 0-50%   |
| <b>CROWN PREP</b>    |                                |                                |                                |         |
| DB1                  | 4-6                            | 4-5                            | ○**                            | 50-100% |
| <b>EXTRACTION</b>    |                                |                                |                                |         |
| ME1, ME2, ME3 ***    | 6                              | 5                              |                                |         |

\* These devices are not provided with a power level indicator. The percentage value indicated in the table refers to the adjustment percentage of the scale of the potentiometer supplied with the device.

\*\* Can be used with or without „soft mode“ function.

\*\*\* Not available in the United States and in Italy market.

# →→ APPROPRIATE INSERT SETTINGS



| function                                | endo | perio / scaler | restorative / prosthetics | power | power |
|---|------|----------------|---------------------------|-------|-------|
| <b>SCALING</b>                          |      |                |                           |       |       |
| S1-S                                    |      | 1-2**          |                           | 1-2   | 1-2   |
| S1, S2, S3, S4, S5, S6                  |      | 1-6**          |                           | 1-6   | 1-4   |
| <b>PERIO</b>                            |      |                |                           |       |       |
| P1, P2, P3, P4, P10, P11, P12, P13, P14 |      | 1-3**          |                           | 1-3   | 1-2   |
| P15, P16R, P16L                         |      | 1-3**          |                           | 1-3   | 1-3   |
| <b>IMPLANT CLEANING</b>                 |      |                |                           |       |       |
| ICS + IC1                               |      | 1-2**          |                           | 1-2   | 1-2   |
| <b>ORTHO CLEAN-UP</b>                   |      |                |                           |       |       |
| D05120, D07102 + D06007                 |      | 3-6            |                           | 3-6   | 2-4   |



| function | endo | perio / scaler | restorative / prosthetics | power | power |
|----------|------|----------------|---------------------------|-------|-------|
|----------|------|----------------|---------------------------|-------|-------|

#### ENDO REVISION

|                    |  |     |  |     |     |
|--------------------|--|-----|--|-----|-----|
| ER1                |  | 1-6 |  | 3-6 | 2-4 |
| ER2, ER3, ER4, ER5 |  | 1-3 |  | 1-3 | 1-2 |

#### ENDO RETRO

|               |  |     |  |     |     |
|---------------|--|-----|--|-----|-----|
| RD3, RD4, RD5 |  | 1-3 |  | 1-4 | 1-3 |
|---------------|--|-----|--|-----|-----|

#### CROWN PREP

|     |  |  |     |     |     |
|-----|--|--|-----|-----|-----|
| DB1 |  |  | 4-6 | 4-6 | 3-4 |
|-----|--|--|-----|-----|-----|

#### EXTRACTION

|                   |  |  |   |  |  |
|-------------------|--|--|---|--|--|
| ME1, ME2, ME3 *** |  |  | 6 |  |  |
|-------------------|--|--|---|--|--|

\* Device no longer available.

\*\* Can be used with or without „soft mode“ function.

\*\*\* Not available in the United States and in Italy market.



# REFERENCE NUMBERS

## → ITEM/REFERENCE NUMBER

|                          |              |
|--------------------------|--------------|
| set maintenance          | 02160035-001 |
| set scaling              | 02160001-001 |
| set perio maintenance    | 02160006-001 |
| set perio treatment      | 02160036-001 |
| set implant cleaning S*  | 02160033-001 |
| set ortho clean-up **    | 02160037-001 |
| set endo revision        | 02160008-001 |
| torque wrench K10        | 02900137-001 |
| dynamometric wrench K6   | 02900074     |
| wrench K1                | 02270003     |
| insert S1-S + wrench K10 | 02960009-001 |
| insert S1 + wrench K10   | 02960001-001 |
| insert S2 + wrench K10   | 02960002-001 |
| insert S3 + wrench K10   | 02960003-001 |
| insert S4 + wrench K10   | 02960004-001 |
| insert S5 + wrench K10   | 02960005-001 |
| insert S6 + wrench K10   | 02960006-001 |
| insert P1 + wrench K10   | 02970001-001 |
| insert P2 + wrench K10   | 02970002-001 |
| insert P3 + wrench K10   | 02970003-001 |
| insert P4 + wrench K10   | 02970004-001 |

## → ITEM/REFERENCE NUMBER



|                               |              |
|-------------------------------|--------------|
| insert P10 + wrench K10       | 03080001-001 |
| insert P11 + wrench K10       | 03080002-001 |
| insert P12 + wrench K10       | 03080003-001 |
| insert P13 + wrench K10       | 03080004-001 |
| insert P14 + wrench K10       | 03080005-001 |
| insert P15 + wrench K10       | 03080006-001 |
| insert P16R + wrench K10      | 03080007-001 |
| insert P16L + wrench K10      | 03080008-001 |
| tip holder ICS                | 03570004-001 |
| tip IC1 (5 pieces)            | 02900112     |
| insert D05120 + wrench K10    | 02990005-001 |
| tipholder D07102 + wrench K10 | 02990007-001 |
| tip D06007 (5 pieces)         | 02900190     |
| insert ER1 + wrench K10       | 03450001-001 |
| insert ER2 + wrench K10       | 03450002-001 |
| insert ER3 + wrench K10       | 03450003-001 |
| insert ER4 + wrench K10       | 03450004-001 |
| insert ER5 + wrench K10       | 03450005-001 |
| insert RD3 + wrench K10       | 03050006-001 |
| insert RD4 + wrench K10       | 03050007-001 |
| insert RD5 + wrench K10       | 03050008-001 |

## → ITEM/REFERENCE NUMBER

|                         |              |
|-------------------------|--------------|
| insert ME1 + wrench K6  | 03130001-001 |
| insert ME2 + wrench K6  | 03130002-001 |
| insert ME3 + wrench K6  | 03130003-001 |
| crown prep set***       | 02160009     |
| tip holder DB1          | 03570001     |
| key AB1                 | 03580001     |
| dynamometric wrench K7  | 02900081     |
| crown prep tip TA12D60  | 03590001     |
| crown prep tip TA14D60  | 03590002     |
| crown prep tip TA16D60  | 03590003     |
| crown prep tip TA12D90  | 03590004     |
| crown prep tip TA14D90  | 03590005     |
| crown prep tip TA16D90  | 03590006     |
| crown prep tip TA14D120 | 03590007     |
| crown prep tip TA16D120 | 03590008     |
| crown prep tip TF12D60  | 03590010     |
| crown prep tip TF16D60  | 03590011     |
| crown prep tip TF12D90  | 03590012     |
| crown prep tip TF16D90  | 03590013     |

\* equipped with: ICS and 5 x IC1 + K10  
 \*\* equipped with: D05120, D07102, 5 X D06007 + 2 K10  
 \*\*\* equipped with: AB1, DB1 + K7

# MECTRON TIP CARD



Ultrasonic inserts wear out the same as manual instruments. Regularly check your inserts by means of the erosion tip card. Whenever your inserts are shorter than the orange hatched line, their performance is 50% less than the initial one. Such inserts should be quickly replaced.



# NOTES

---



mectron s.p.a., via Loreto 15/A,  
16042 Carasco (Ge), Italy,  
tel +39 0185 35361, fax +39 0185 351374

→ [www.mectron.com](http://www.mectron.com) – [mectron@mectron.com](mailto:mectron@mectron.com)



© Copyright mectron S.p.A., Carasco, Italy  
All rights reserved. Text, images and graphics contained in Mectron brochures are protected by copyright and other property rights laws. Content cannot be copied, disclosed, modified or made accessible to third parties for commercial purposes without written permission from Mectron S.p.A.

The images shown are for illustration purposes only and may not be an exact representation of the product. Please consult product labels and manual for indications, contraindications, hazards, warnings, cautions. Please contact your local sales representative for product availability.

The following products are registered under MDR 2017/745. For graphic reasons the styling of the product names might be different. The registered product names correspond as following:

|               |                    |                     |                   |
|---------------|--------------------|---------------------|-------------------|
| • combi touch | COMBI <i>touch</i> | • COMPACT PIEZO LED | COMPACT PIEZO LED |
| • multipiezo  | MULTIPIEZO         | • COMPACT PIEZO P2K | COMPACT PIEZO P2K |
| • MICROPIEZO  | MICROPIEZO         |                     |                   |