

**DENTSPLY**  
MAILLEFER

CHEMIN DU VERGER 3  
CH-1338 BALLAIGUES  
SWITZERLAND

[www.dentsplymaillefer.com](http://www.dentsplymaillefer.com)

CE  
0459

INL/F19 03 54.X/00/2013 - created 01/2013

**DENTSPLY**  
MAILLEFER

+  
WE  
KNOW  
ENDO.



**PROTAPER • NEXT™**  
flexible performance



### GLIDE PATH



P1

PATHFILE®



P2

### SHAPING



X1



X2

PROTAPER•NEXT™

### OPTIONAL



X3



X4



X5

## PATHFILE® GLIDE PATH

P1



013 - 2%

P2



016 - 2%

## PROTAPER•NEXT™ SHAPING

X1



017 - 4%\*

X2



025 - 6%\*

X3



030 - 7%\*

X4



040 - 6%\*

X5



050 - 6%\*

\*variable taper