

X-Smart® Apex Locator

The right way to the apex





X-Clusive technology for X-Treme precision

X-Smart® Apex Locator offers highly reliable accuracy—with patented measurement technology, an easy-to-use touchscreen interface, and audio/visual cues for added guidance during procedures.





Measuring with X-Pert accuracy

Other apex locators use mixed AC signal frequencies to detect the apex; X-Smart® Apex Locator alternates them, canceling the need for signal filtering and eliminating the noise that can interfere with performance.



Features designed for you

- 1 Large, full-color, high-resolution touchscreen that works with gloves
- 2 Dr's Choice feature allows for saved user preferences and provides visual and audio cues when the file tip has reached the pre-selected position
- 3 Visual and audio guidance during measurements
- 4 Long-lasting, lithium-ion rechargeable battery with fast charging
- 5 USB Type-C connector for more durability and flexible orientation
- 6 Connection of File Clip and Lip Clip to any socket of the measurement cable with onboard back storage

Get up to **150** root canal treatments with a full charge*



An accurate working length determination is crucial to the success of your root canal treatment.

Let X-Smart® Apex Locator guide your way. Learn more at dentsplysirona.com/XSmartApexLocator

Technical Specifications

- 5" touchscreen interface
- Lithium-ion battery
- Fully recharged in approximately 3 hrs
- On/Off button
- Sound control
- Device settings
- Automatic turn-off after 5-15 min
- Dimensions: Width 95 mm x Depth 146 mm x Height 111 mm
- Calibration not required
- Frequencies used for measurements 500 Hz and 8 kHz
- Product Warranty : 24 months / Accessories Warranty : 6 months



Product Reference

X-Smart® Apex Locator		ACCESSORIES
	B00XSAL000000	1x X-Smart® Apex Locator
		1x Measurement Cable
		2x File Clip
		2x Lip Clip
		1x Touch Probe
		1x Charger with Plugs (EUR, UK, US, AUS)
		B00XSALAC1000 X-Smart® Apex Locator Measurement Cable
		B00XSALAC2000 X-Smart® Apex Locator File Clip, set of 2
		B00XSALAC3000 X-Smart® Apex Locator Lip Clip, set of 5
		B00XSALAC4000 X-Smart® Apex Locator Touch Probe, set of 2
		B00XSALAC5000 X-Smart® Apex Locator Charger with Plugs (EUR, UK, US, AUS)