

# MULTIMEDIA TROLLEY




## WHY MULTIMEDIA TROLLEY

- All-in-one working station (Iris View light, medical monitor, recording and streaming device)
- Compact dimensions, on wheels
- Power supply 230V
- Wi-Fi and Ethernet for internet and local Network connection

## WHY G.COMM



We focus on R&D  
and we apply innovative  
technologies



We are proud to be a manufacturer  
and our whole production  
chain is in Italy



We are close to customers  
and support them with  
qualified technical assistance




[www.gcomm-online.com](http://www.gcomm-online.com)

- Shows your video directly on monitors, smartphones, tablets or PCs.
- Records your videos on internal storage, a connected USB drive or an associated smartphone or tablet.
- Live streaming to popular services such as Facebook Live, YouTube, Twitch or to a RTMP/ RMTSP server .
- Internet connectivity via wired Ethernet, Wi-Fi and optional 3G/4G modem support.



LIVE SURGERY  
STREAMING EASY FOR  
EVERYONE



**Specific App to see  
your videos on iOS  
or Android devices**

**Record videos or  
take pictures directly  
on a USB Drive or  
smartphones**

**It creates a video  
stream to your local  
PCs or a LIVE  
stream to Youtube e  
facebook**

